

Analysis Date: 13-05-2025

Re-test Date: 13-05-2028

## AFATINIB SULFONIC ACID ANALOG

### Identification


**Chemical Name** : 4((4-((3-Chloro-4-fluorophenyl) amino)-7-(((S)-tetrahydrofuran-3-yl) oxy) quinazolin-6-yl)amino)-1-(dimethylamino)-4-oxobutane-2-sulfonic acid

**CAT No** : ALL-A08000

**CAS No** : 2283295-16-7

**Molecular Formula** : C<sub>24</sub>H<sub>27</sub>ClFN<sub>5</sub>O<sub>6</sub>S

**Molecular Weight** : 568.0



### Analytical Information

|                              |                   |                          |           |
|------------------------------|-------------------|--------------------------|-----------|
| <b>Batch Code</b>            | : ALL-A08000      | <b>HPLC Purity</b>       | : 98.02%  |
| <b>Solubility</b>            | : MeOH: ACN       | <b>Potency</b>           | : 95.80%  |
| <b>Appearance of Product</b> | : Off White Solid | <b>Mass</b>              | : Confirm |
| <b>Long Term Storage</b>     | : 2-8 °C          | <b>IR Analysis</b>       | : Confirm |
| <b>Weight Content By TGA</b> | : 1.255%          | <b><sup>1</sup>H NMR</b> | : Confirm |
| <b>Residue Of Ignition</b>   | : 1.009%          |                          |           |

### Additional Information

%Potency = [100 - (Inorganic Impurities % + Water %) × Chromatographic Purity%]/100 =

[100 - (1.255+1.009 × 98.03)]/100 = 95.80%

**Recommendation** : Released

|                          | Department | Name                               | Signature |
|--------------------------|------------|------------------------------------|-----------|
| Prepared and Reviewed by | Analytical | Mr. Vipul Khadse<br>Jr. Executive  |           |
| Approved By              | QA&QC      | Dr. Ashish Keche<br>Director QA&QC |           |

**Attachment** : HPLC, Mass, <sup>1</sup>H NMR, IR, TGA

**Shipping Condition** : All Product are stable to be shipped at room temperature, unless otherwise specified

### Corporate Office

Allmpus Laboratories Pvt Ltd 38, Solitaire Industrial Estate, Plot No.TS-1, Near Nandi Palace Hotel, MIDC-Phase 1, Dombivli (E), Mumbai, Maharashtra 421203