

Analysis Date:

Re-test Date:

TESTOSTERONE-D8

IDENTIFICATION


Chemical Name : (8R,9S,10R,13S,14S,17S)-17-hydroxy-10,13-dimethyl-1,2,6,7,8,9,10,11,12,13,14,15,16,17-tetradecahydro-3H-cyclopenta[a]phenanthren-3-one-2,2,4,6,6,16,16,17-d8

CAT No : ALL-D0306

CAS No : 2416308-73-9

Molecular Formula : C₁₉H₂₀D₈O₂

Molecular Weight : 296.5



Analytical Information

Batch Code : ALL-D0306

Solubility : MeOH: ACN

Appearance of Product : Off White Solid

Long Term Storage : 2-8°C

Weight Loss By TGA : N.A

Residue Of Ignition : N.A

HPLC Purity : NLT 98.00 %

Potency : NLT 97.00 %

Isotopic Enrichment : NLT 98.00 %

IR Analysis : Confirm

Mass : Confirm

¹H NMR : Confirm

Additional Information

%Potency = [100 - (Weight Loss By TGA % + Residue Of Ignition %) × Chromatographic Purity%]/100 = N.A

Recommendation : Released

	Department	Name	Signature
Prepared and Reviewed by	Analytical	Mr. Vipul Khadse QA Officer	
Approved By	QA&QC	Dr. Ashish Keche Director QA&QC	

Attachment : HPLC, Mass, ¹H NMR, IR, TGA

Shipping Condition : All Product are stable to be shipped at room temperature otherwise specified

Corporate Office

Allmpus Laboratories Pvt Ltd 38, Solitaire Industrial Estate, Plot No.TS-1, Near Nandi Palace Hotel, MIDC-Phase 1, Dombivli (E), Mumbai, Maharashtra 421203